

FREESTYLE - Grawery na soczewkach

| Nazwa soczewki | Index | Grawer od strony nosa |
|---|-------|-----------------------|
| Freestyle People 1.50 CR-39 | 1,50 | PLC |
| Freestyle People 1.53 TVX | 1,53 | PLV |
| Freestyle People 1.60 MR-8 | 1,60 | PLM |
| Freestyle People 1.67 MR-6 | 1,67 | PLH |
| Freestyle People 1.50 Transitions | 1,50 | PLA |
| Freestyle People 1.60 Transitions | 1,60 | PLB |
| Freestyle People 1.67 Transitions | 1,67 | PLB |
| Freestyle People 1.53 Transitions TVX | 1,53 | PFV |
| Freestyle People 1.50 Polarized | 1,50 | PLD |
| Freestyle People 1.60 Polarized | 1,60 | PLE |
| Freestyle People 1.74 | 1,74 | PLU |
| Freestyle Office 1.50 CR-39 | 1,50 | OFC |
| Freestyle Office 1.53 TVX | 1,53 | OFV |
| Freestyle Office 1.60 MR-8 | 1,60 | OFM |
| Freestyle Office 1.67 MR-6 | 1,67 | OFH |
| Freestyle Office 1.50 Transitions | 1,50 | OFA |
| Freestyle Light 1.50 CR-39 | 1,50 | FLC |
| Freestyle Light 1.60 | 1,60 | FLM |
| Freestyle Light 1.67 | 1,67 | FLH |
| Freestyle Light 1.56 Sunsenors | 1,56 | FLS |
| Freestyle Intelligent 1.50 CR-39 | 1,50 | DIC |
| Freestyle Intelligent 1.53 TVX | 1,53 | DIV |
| Freestyle Intelligent 1.60 MR-8 | 1,60 | DIM |
| Freestyle Intelligent 1.67 MR-6 | 1,67 | DIH |
| Freestyle Intelligent 1.50 Transitions | 1,50 | DIA |
| Freestyle Intelligent 1.60 Transitions | 1,60 | DIB |
| Freestyle Intelligent 1.67 Transitions | 1,67 | DIT |
| Freestyle Intelligent 1.53 Transitions TVX | 1,53 | DFV |
| Freestyle Intelligent 1.53Nxt Sun/Varia/Polar | 1,53 | DIN |
| Freestyle Intelligent 1.50 Polarized | 1,50 | DID |
| Freestyle Intelligent 1.60 Polarized | 1,60 | DIE |
| Freestyle Intelligent 1.74 | 1,74 | DIU |